



UltraFast
Innovations

YOUR KEY to innovation and success



XUV Phase Retarder **AURORA**


Our XUV Phase Retarder acts as a quarter waveplate to turn linearly polarized XUV light into circularly polarized light without introducing noticeable dispersion. The phase retarder achieves close-to-circular polarization of $P_c=0.75$ and obtains $> 25\%$ transmission around 66 eV photon energy, where the Ni M_2/M_3 edge locates [1].

A broad spectral range from 40 to 85 eV is supported to cover the M_2/M_3 edge of 3 transition metals - Fe, Co and Ni.

The retarder uses a transmission optimized, four mirror grazing incidence reflection geometry that induces a quarter wave phase offset between the s- and p-polarization components of a linearly polarized input XUV beam. At the Ni M_2/M_3 transition up to 3% dichroism contrast is observed (i.e. $> 85\%$ of the theoretical value) [1].

A clear aperture of 3 mm will allow the low divergent XUV to pass through without clipping. Our XUV

Phase Retarder is ideally geared to be combined with ultrafast high-harmonic XUV sources adding spin-sensitivity to conventional laser based pump-probe experiments via attosecond magnetic circular dichroism detection.



Key Product Features:

- Create circularly polarized XUV light without adding dispersion (ideal for attosecond applications)
- Up to $\approx 40\%$ max transmission
- Broad spectral range covering 40 – 85 eV
- Compact design, with minimum change in the delay line
- Entrance and exit iris diaphragms for alignment
- Easy to implement in existing setups

Characteristics:


Spectral Range	40-85 eV
Avg Transmission	$>25\%$
Max Ellipticity P_c^1	0.75 @ 66 eV
Extra Beam Path	3 mm
Clear Aperture	3 mm
Footprint	40mm (D) x 92mm (L)

¹Degree of Ellipticity expressed as Stokes Parameters $P_c=S3/S0$


UltraFast Innovations GmbH
Dieselstr. 5
85748 Garching
Germany

phone: +49 89 36039 - 437
fax: +49 89 36039 - 453
info@ultrafast-innovations.com
www.ultrafast-innovations.com





Phase and Transmission, p-Pol



Phase and Transmission, s-Pol

Element	Fe	Co		Ni	
Edge	M_2/M_3	M_2	M_3	M_2	M_3
Energy	52.7 eV	58.9 eV	59.9 eV	66.2 eV	68 eV
Transmission	9%	17%	18%	26%	30%

References:

- [1] Siegrist, F., Gessner, J.A., Ossiander, M. et al. Light-wave dynamic control of magnetism. *Nature* **571**, 240-244(2019).